

EBM Mouth Care

Scenario

Mouth care is being prepared using Orders. You have logged into the back office.

OR

Mouth care is being administered. You have logged into the handheld device.

Instruction: Back Office

Mapp | Mobile App Platform

PharmaTrac
Safe medication administration

MilkTrac
Safe feeding of babies

Dashboard
Baby associations
Batch correction
Discharge
Orders
Batches
Home labels
Reports
Management

MealTrac
Meal management

Repacker
Undose preparation

Dispenser
Medication dispensing

Configuration

Orders

Search:

New preparation

Process	Baby	Birth date	Ward	Bed	Prescriptions	Milk type	Collected (ml)	Stored (ml)	Prepared (ml)
2012854	Don Norma Beth	2021-06-07	Medicine - East Wing	228	35 ml 3/3h	Maternal Milk		200	
9199959	DUMMY PATIENT2 2012854 Doris	2021-06-07	Medicine East Wing					200	
da	DUMMY PATIENT da	2019-12-02	Paediatric						
db	DUMMY PATIENT db	2019-12-02	Paediatric						
dc	DUMMY PATIENT dc	2019-12-02	Paediatric						
dd	DUMMY PATIENT dd	2019-12-02	Paediatric						
1030858	DUMMY PATIENT 1030858 TERESA ARRUDA BOTELHO	2021-08-21	Intensivists - Neonatology						
11223344	Dummy Teste Labrator Dashboard	2022-07-28	Medicine East Wing						
123456	DUMMY PATIENT 123456 SURNAME	2000-01-01	Nephrology	1					
123456	DUMMY PATIENT 123456 SURNAME2	2000-05-05	Nephrology	2					
2289002	DUMMY PATIENT 2289002 Williams	2021-06-01	Medicine - East Wing	220	15 ml 3/3h	Maternal Milk			
2890866	DUMMY PATIENT 2890866 Diaz	2021-06-04	Medicine - East Wing	220	40 ml 2/2h	Maternal Milk			
3006313	Mary Stone	1945-02-02	Medicine - East Wing	8	Unrestricted	Normal			
5006605	Sam Brutus Afrikaans	1940-10-03	Medicine - East Wing	A		Pasty			
751084	DUMMY PATIENT 751084 BIG	2021-05-01	BIG EtherOx						

Showing from 1 to 15 Total: 25 records

Previous 1 2 Next

1. Select **Orders**
2. Search for your patient – you can search by MRN
3. Select the patient
4. Click **New preparation**
5. The **New preparation** window will now open

Platform

New preparation

BABY NAME
DUMMY PATIENT2 2012854 Doris

MILK TYPE

Add **Clear list** **Validate**

Type	Exp. date	Batch	Quantity (ml)	Prep. date	Stored (ml)	Prepared (ml)
Maternal Milk	2023-03-05 14:00	L380	1			

Automatic segmentation

Preparations Num. 1 Quantity (ml) 1.00

Confirm

Administration history

Date	Type	Administered (ml)	Remaining (ml)
2022-09-14 12:34	Maternal Milk	50	0
2022-09-14 12:35	Maternal Milk	200	0
2022-09-08 15:17	Maternal Milk	50	0

Stock


Expiry date	Quantity (ml)	Batch	Status	fridge
2023-03-05 14:00	200	L380	Frozen	Practically

6. Click **Add**

7. The **Add milk** window should now be open
8. Choose the desired batch (the batches are sorted by expiration date)
9. Choose the checkbox next to **Mouth care**. This defaults to 1ml, but you can adjust the quantity
10. Select **Add milk**
11. You should then be returned to the **New preparation** window

12. Click **Validate**
13. The **Validate** window will open:
14. Two labels will be printed, as two new batches have been created. One label should be attached to the mouth care, and the second label should be attached to the remaining EBM which should then be stored.
15. Click **Finish**

Instruction: Handheld Device



1. Select the **Options Menu** on the handheld device (this is the fourth screen)
2. Select **Mouth Care**

3. The **Mouth Care** screen should appear as above:
4. Scan the EBM
5. Scan the patient's wristband
6. If there is a mismatch for either the mouth care, or the wristband, an error message will appear.